Quantitation and Ratio Determination of Uranium Isotopes in Water and Soil Using Inductively Coupled Plasma Mass Spectrometry (ICP-MS)

D.N. Kurk, T.E. Beegle, S.C. Spence and R.J. Swatski



maintaining the data needed, and c including suggestions for reducing	lection of information is estimated to completing and reviewing the collect this burden, to Washington Headqu uld be aware that notwithstanding ar DMB control number.	ion of information. Send comments arters Services, Directorate for Infor	regarding this burden estimate or mation Operations and Reports	or any other aspect of the 1215 Jefferson Davis	is collection of information, Highway, Suite 1204, Arlington	
1. REPORT DATE MAR 2012		2. REPORT TYPE		3. DATES COVE 00-00-2012	red 2 to 00-00-2012	
4. TITLE AND SUBTITLE		5a. CONTRACT NUMBER				
Quantitation and Ratio Determination of Uranium Isotopes in Water and Soil Using Inductively Coupled Plasma Mass Spectrometry (ICP-MS)					5b. GRANT NUMBER	
Son Osing muuch	refy Coupled I fasilia	5c. PROGRAM ELEMENT NUMBER				
6. AUTHOR(S)		5d. PROJECT NUMBER				
					5e. TASK NUMBER	
				5f. WORK UNIT NUMBER		
	ZATION NAME(S) AND AD Public Health,5158 I ID,21010-5403		berdeen	8. PERFORMING REPORT NUMB	GORGANIZATION ER	
9. SPONSORING/MONITORING AGENCY NAME(S) AND ADDRESS(ES)				10. SPONSOR/MONITOR'S ACRONYM(S)		
				11. SPONSOR/M NUMBER(S)	ONITOR'S REPORT	
12. DISTRIBUTION/AVAIL Approved for publ	LABILITY STATEMENT ic release; distributi	on unlimited				
13. SUPPLEMENTARY NO Presented at the 9t 26-29 March 2012	h Annual DoD Envi	ronmental Monitori	ing and Data Qua	ality (EDMQ)) Workshop Held	
14. ABSTRACT						
15. SUBJECT TERMS						
16. SECURITY CLASSIFIC	17. LIMITATION OF	18. NUMBER	19a. NAME OF			
a. REPORT unclassified	b. ABSTRACT unclassified	c. THIS PAGE unclassified	Same as Report (SAR)	OF PAGES 17	RESPONSIBLE PERSON	

Report Documentation Page

Form Approved OMB No. 0704-0188





Disclaimer

The views expressed herein are the views of the author(s) and do not reflect the official policy of the Department of the Army or the Department of Defense, or the U.S. Government.

Use of trademarked name does not imply endorsement by the U.S. Army but is intended only to assist in identification of a specific product.





Introduction

- Uranium overview
- **♦ Sample prep:**

water – EPA 3020

Soils – EPA 3052 (modified)

♦ Analysis – ICP-MS

water - EPA 200.8

soils - **EPA** 6020

Quality Controls

- ♦ ICP-MS and α-Spec
- **♦** Summary
- Questions





Introduction

- Ubiquitous element
- Naturally Occurring U Isotopes:
 234 (0.0055%, 0.245 E6 yr.), 235 (0.72%, 703 E6 yr.), 238 (99.275%, 4,468 E6 yr.)
- **♦ Natural U235/238 <u>atomic</u> Ratio: 7.2 x 10-3**
- ♦ Natural U234/238 α activity Ratio: 1 (secular equilibrium)
- **♦ Used for fuel in atomic energy and warfare**
- → Depleted Uranium DU: 235 Isotope Quantity Reduced U235/238 atomic Ratio: 2 x 10-3





Sample Preparation

- Water EPA 3020
 Acid digestion
- ♣ Soils EPA 3052 (modified)
 Acidic microwave digestion
 Complete digestion
- → Ratio Depends on Matrix (see methods above)





Sample Analysis ICP-MS

- **♦ Water EPA 6020**
- **♦ Soils EPA 200.8**
- → Ratios In house Method





Sample Analysis QC

- **♦ Sample Duplicates precision check**
- **♦ Blanks contamination check**
- **♦ Laboratory Control Samples accuracy check**
- **♦ Matrix Spikes matrix effect**
- *Mass bias correction standard





Analysis Recovery QC

	Water	Soil	Ratio
Duplicates		≤ 20 % RPD	
Blanks	< RL	< RL	
LCS	± 15 %	± 20 %	
MS	± 30 %	± 30 %	
Inst Spike	± 20 %	± 20 %	
ISA / ISB		± 20 %	





Common Analysis Techniques

- **♦** α Spectroscopy
- **♦ ICP-MS**





a Spectroscopy

- Measures 234 and 238 isotopes
 U-234 from the Uranium Decay Series
 U238 → Th234 → Pa234 → U234 → Th230 →
- **♦** Sample preparation required (matrix removed)
- **♦** Tracer added for quantification
- **Measure** α particles from radioactive decay
- ♦ Ratio and Concentration in same analysis
- ♦ Detection limits depends on count time





ICP-MS

- Measures 235 and 238 isotope ions
 235 from the Actinium Decay Series
 U235 → Th231 → Pa231 → Ac227 →
- **♦** Sample preparation
- **♦** Count ions
- **♦ Conc. and Ratio: two different analyses**
- **♦** Detection limit: matrix and instrument





ICP- MS 235/238 Ratio Comparison

Uncorrected Bias	Corrected Bias
6.51 x10-3	7.24 x10-3
6.82 x10-3	7.18 x10-3
6.58 x10-3	7.25 x10-3
6.24 x10-3	7.14 x10-3
6.74 x10-3	7.22 x10-3

Accepted Ratio value 7.26 x10-3





ICP-MS

- Measures 235 and 238 isotope ions
 235 from the Actinium Decay Series
 U235 → Th231 → Pa231 → Ac227 →
- **♦ Sample preparation**
- **♦** Count ions
- **♦ Conc. and Ratio: two different analyses**
- **♦** Detection limit: matrix and instrument





Choices

- **♦** α Spectroscopy
- **♦ ICP-MS**

♣ Questions you need answered: Concentration? Ratio?





ICP- MS and α Spec

ICP-MS	α Spec
235 Conc. Sufficient for ratio	Long count time
Ratio & Conc. Separate Analysis	Ratio & Conc. Same Analysis
Correct ratio? – bias, conc.	Correct ratio – recoil effect
100 mL sample	1 L sample





Acknowledgments

Army Institute of Public Health (AIPH) EMDQ





Questions

???